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### PRODUCT DESCRIPTION

**Intended use**

The NycoCard<sup>®</sup> HbA1c Control contains liquid preparations of stabiliseried bovine haemoglobin.

Quality control using the NycoCard<sup>®</sup> HbA1c should be performed to confirm that the NycoCard<sup>®</sup> HbA1c test kit and NycoCard<sup>®</sup> READER II are working properly. Only when controls are used routinely, and the values obtained are within the acceptable ranges, can accurate results for patient samples be assured.

**Kit contents**

1 Normal Control (1.5 mL). Stabilised bovine haemoglobin preparation.

1 Abnormal Control (1.5 mL). Stabilised bovine haemoglobin preparation.

1 Package insert

**Warnings and precautions**

- For *in vitro* diagnostic use.
- Use gloves when handling NycoCard<sup>®</sup> HbA1c control solutions.
- No known test method can offer complete assurance that products derived from animal blood will not transmit disease. Such materials should always be handled as potentially infectious.
- Do not use the NycoCard<sup>®</sup> HbA1c controls after their expiry date or if they have not been stored in accordance with recommendations.
- Discard the vial if there is evidence of microbial or fungal contamination.

### STABILITY AND STORAGE

- Avoid exposure to direct sunlight and temperatures above 8°C.
- Do not freeze.

**Unopened control vials**

The expiry date of the control kit only applies if the product is stored refrigerated (2-8°C) in the original container. The expiry date is the last day of the month stated on the vial label.

**Opened control vials**

- Opened control vials are stable for thirty (30) days when stored refrigerated. Note the date of opening on the vial label.
- Always store the control vial refrigerated when not in use.
- Opened control vials should be stored in an upright position.

### TEST PROCEDURE

**Important procedural notes!**

- Always store the control vials in a refrigerator (2-8°C) in an upright position.
- Clean the outside of the control vial neck and replace the cap after use.
- Immediately return the control back to the refrigerator.

**Procedure**

- Remove NycoCard<sup>®</sup> HbA1c Control vial(s) from refrigerator and equilibrate to room temperature for a minimum of 30 minutes.
- Shake vial vigorously for a minimum of ten (10) seconds until a uniform suspension is obtained.
- Remove cap from vial and then tilt vial slowly until suspension is in neck of vial.
- Fill the capillary by touching the control suspension with the end of the capillary. Ensure that the capillary is completely filled.
- Clean the outside of the control vial neck and replace the cap.

The NycoCard<sup>®</sup> HbA1c Control should be assayed in the same manner as ordinary patient samples. Follow the procedure described in the package insert, enclosed in each NycoCard<sup>®</sup> HbA1c test kit.

### EXPECTED VALUES

NycoCard<sup>®</sup> HbA1c is traceable to the IFCC Reference Method<sup>1</sup>. The values reported are DCCT/NGSP aligned.

An acceptable range for analysis using NycoCard<sup>®</sup> HbA1c and NycoCard<sup>®</sup> READER II has been assigned for these control materials. The measured value should be within this acceptable range. Re-test the control material if the control result is outside the range. Re-test a sample from a new control vial if the measured value is still outside the acceptable range. Patients' test results must be declared invalid if the control does not perform as expected. If the problem persists, see "Troubleshooting" in the NycoCard<sup>®</sup> HbA1c package insert or contact your local supplier.

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