NycoCard[®] **HbA1c**

REF

1042184 : 24 Tests

CE

GB	3 ... 6
DE	7 ... 10
FR	11 ... 14
ES	15 ... 18
IT	19 ... 22
NL	23 ... 26
GR	27 ... 30
NO	31 ... 34
SE	35 ... 38
DK	39 ... 42
FI	43 ... 46

AVAILABLE PRODUCTS

REF	1042184 NycoCard® HbA1c	24 tests
REF	1113022 NycoCard® HbA1c Control	2 x 0.4 mL
REF	1113117 NycoCard® READER II	1 instrument
REF	1054331 Capillary tubes, 5 µL	200 units
REF	1003668 Capillary tube holder	1 unit
REF	1113060 Test tube rack (for 6 tubes)	1 unit
REF	1113059 Mini-Pet stand	1 unit
REF	1017483 Mini-Pet, 25 µL	1 unit
REF	1040711 Pipette tips, 2-200 µL	2 x 500 units

PRODUCT DESCRIPTION

Intended use

NycoCard[®] HbA1c is a rapid *in vitro* test for the measurement of glycated haemoglobin in human blood.

Test principle

NycoCard[®] HbA1c is a boronate affinity assay. The kit contains test devices with a porous membrane filter, test tubes pre-filled with reagent and a washing solution. The reagent contains agents that lyse erythrocytes and precipitate haemoglobin specifically, as well as a blue boronic acid conjugate that binds cis-diols of glycated haemoglobin.

When blood is added to the reagent, the erythrocytes immediately lyse. All haemoglobin precipitates. The boronic acid conjugate binds to the cis-diol configuration of glycated haemoglobin.

An aliquot of the reaction mixture is added to the test device, and all the precipitated haemoglobin, conjugate-bound and unbound, remains on top of the filter. Any excess of coloured conjugate is removed with the washing solution.

The precipitate is evaluated by measuring the blue (glycated haemoglobin) and the red (total haemoglobin) colour intensity with the NycoCard[®] READER II, the ratio between them being proportional to the percentage of HbA1c in the sample.

Kit contents, 24 test kit

TD/Test Device 1 x 24 units

Plastic device containing a membrane filter.

R1/Reagent 1 x 24 x 0.2 mL

Glycinamide buffer containing Zn ions, dye-bound boronic acid and detergents.

R2/Washing Solution 1 x 2.0 mL

Morpholine buffered NaCl solution and detergents.

Materials required

(not supplied with the kit)

- Capillary tube or pipette (5 µL) for sampling.
- Pipette (25 µL) and pipette tips for the application of reaction mixture and R2/Washing Solution.
- NycoCard[®] READER II for measurement of the test result.

Warnings and precautions

IVD For *in vitro* diagnostic use only.

- R1/Reagent and R2/Washing Solution contain sodium azide (0.05%), which is a toxic agent.
- R1/Reagent contains formamide (6.2%), which is irritating to eyes and skin, and is a toxic agent.
- R2/Washing Solution contains morpholine (50 mmol/L), which is a toxic agent.

Analytical specificity

NycoCard[®] HbA1c measures the total glycated haemoglobin (GHb), but reports a standardised HbA1c value.

Standardisation

Axis-Shield PoC AS works in close relationship with the European Reference Laboratory for Glycohemoglobin (ERL). Standardisation of NycoCard[®] HbA1c at DCCT level is carried out according to IFCC recommendations.

Measuring range

Measuring range: 3-18% HbA1c

Measuring interval: 0.1% HbA1c

Reference range for NycoCard[®] HbA1c

4.5-6.3% HbA1c

Precision

The coefficient of variation (CV) is below 5% in the whole measuring range.

Limitations of the test

- Elevated amounts of glucose, bilirubin, lipids and fructosamine were added to blood samples with normal (5.1%) and

elevated (9.7%) HbA1c values. No interference was obtained.

- Pre-glycated haemoglobin does not interfere with the test.
- Haemolysed samples with plasma Hb >3 g/100 mL will interfere with the test system.
- NycoCard® READER II corrects for Hb-concentrations in the range 6-18 g/100 mL.

STABILITY AND STORAGE

Unopened kits

The expiry date of the kit applies to storage at 2-8°C in original container. Exposure of kit components to temperatures above 25°C or humidity above 70% should be avoided. Do not freeze.

Opened kits

- **R1/Reagent** must be stored dark at 2-8°C. Equilibrate the R1/Reagent to room temperature (20-25°C) before use. Equilibration to room temperature can be achieved by holding the tube in a closed hand for 30 seconds. The R1/Reagent tube can be stored for maximum 6 hours at room temperature before use. Avoid direct sunlight.
- **TD/Test Device** can be stored at room temperature (15-25°C). Store the TD/Test Devices in the original bag and avoid humidity below 20% and above 70%. The TD/Test Device should have room temperature when used.
- **R2/Washing Solution** can be stored at room temperature (15-25°C). The R2/Washing Solution should have room temperature when used.

Sample material

Blood samples can be stored up to 10 days at 2-8°C before analysis. Avoid measuring haemolysed samples (see paragraph "Limitations of the test"). Do not freeze.

TEST PROCEDURE

Important procedural notes!

- Do not interchange components from different kits or kit lots.
- Equilibrate the R1/Reagent to room temperature (20-25°C) before use.
- Do not touch the membrane with the pipette tip.
- Change the pipette tip between each pipetting step.

Sample material

Capillary blood and venous blood with or without anticoagulant (EDTA, heparin and NaF) can be used.

Internal quality control

A quality control material with NycoCard® HbA1c specific target values should be used to confirm the efficacy of the reagents and the correct performance of the test. Lyophilised control materials should not be used with the NycoCard® HbA1c assay.

1 Precipitation of haemoglobin

Add 5 µL whole blood to the test tube pre-filled with R1/Reagent. Mix well. Leave the tube for minimum 2 minutes, maximum 3 minutes.

Note! Equilibrate the R1/Reagent to room temperature (20-25°C) before use.

2 Application of sample

Remix to obtain a homogenous suspension. Apply 25 µL of the reaction mixture to a TD/Test Device by holding the pipette approx. 0.5 cm above the test well. Empty the pipette quickly in the middle of the test well. Allow the reaction mixture to soak completely into the membrane (approx. 10 seconds).

Note! Avoid air bubbles.

3 Application of R2/Washing Solution

Apply 25 µL R2/Washing Solution to the TD/Test Device. Allow the washing solution to soak completely into the membrane. Wait for minimum 10 seconds.

Note! Avoid air bubbles.

4 Test result measurement

Read the test result within 5 minutes using NycoCard® READER II.

Note! Further instructions are given in the NycoCard® READER II instruction manual.

TROUBLE SHOOTING

Problem	Possible cause	Corrective action
The precipitate does not cover the whole surface of the test hole (white spots on the membrane surface).	The reaction mixture has been added too slowly or on the "wall" of the test well.	Retest the sample, and add the reaction mixture quickly in the middle of the test well.
	An air bubble has been formed when applying the reaction mixture to the test well.	Retest the sample and avoid air bubbles.
	Incorrect storage of TD/Test Device.	Retest using correctly stored TD/Test Device.
No precipitate on the TD/Test Device membrane.	Incomplete mixing before incubation, i.e., the capillary tube was not empty.	Retest the sample. Mix the tube more thoroughly.
The NycoCard® READER II displays "Hb conc. too low".	The haemoglobin concentration of the sample is below the measuring range.	Retest the sample and use 10 µL of whole blood: fill 2 capillary tubes simultaneously and drop both capillaries into one test tube with R1/Reagent. The NycoCard® READER II will display the correct test result.
The NycoCard® READER II displays "Reduce Hb conc.".	The haemoglobin concentration of the sample is above the measuring range.	Mix the R1/Reagent from two test tubes (400 µL) and retest the sample (5 µL). The NycoCard® READER II will display the correct test result.
	The volume of R2/Washing Solution added to the test device is too small.	Retest the sample. Make sure that a correct volume of R2/Washing Solution is applied (25 µL).
Unexpected low results.	The R1/Reagent has not been equilibrated to room temperature (20-25°C) before use.	Retest the sample. Make sure that the R1/Reagent is equilibrated to room temperature before use.
	The R1/Reagent has been stored at room temperature/not protected from sunlight.	Retest the sample with a R1/Reagent stored refrigerated (2-8°C)/protected from sunlight.
	The sample material used was haemolysed.	Retest using fresh sample material.
The R2/Washing Solution does not soak into the TD/Test Device membrane.	The R1/Reagent has not been equilibrated to room temperature (20-25°C) before use.	Retest the sample. Make sure that the R1/Reagent is equilibrated to room temperature before use.
	The sample material used was haemolysed.	Retest using fresh sample material.

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Medinor Finland Oy Ab
Axis-Shield Diagnostics Ltd.