ALPHADOTS GANGLIOSIDES PROFILE

PROFILE IgG and/or IgM
Code : AD GANG12 – 12 tests

PROFILE IgG + IgM
Code : AD GANGS12 – 12 tests
AD GANGS24 – 24 tests

DEUTSCH (D)

ENGLISH (GB)

ESPAÑOL (ES)

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ENGLISH

IMMUNODOT ASSAY FOR THE DETERMINATION OF hIgG & hIgM OR hIg (G+M) AUTOANTIBODIES AGAINST GANGLIOSIDES IN HUMAN SERUM.

AD-GANG12, AD-GANGS12 12 Determinations
AD-GANGS24 24 Determinations

FOR IN VITRO DIAGNOSTIC USE ONLY

1. CLINICAL APPLICATIONS

Gangliosides

Gangliosides are sialic acid containing glycopingolipids composed of a long-chain aliphatic amine (ceramide) attached between one and five hexoses, at least one of which must be of sialylated. The presence of sialic acid molecule(s) attached to galactose residue(s) in the hexose core defines as glycopingolipid as a ganglioside (1).

Gangliosides are localized in the outer layer of plasma membranes and are abundant in the myelin sheath of Schwann cells (peripheral nervous system) and of oligodendrocytes (central nervous system) (2).

Neuropathy and anti-Glycoconjugate autoantibodies

The first clinical syndrome, in which anti-Glycoconjugate autoantibodies were described was the IgM Paraproteinemic Neuropathy associated with anti-MAG autoantibodies (1). Further studies revealed the presence of Peripheral Neuropathy and anti-Glycolipid autoantibodies such as anti-Ganglioside or anti-MAG in a large portion of peripheral neuropathic conditions (for review: 2-6). Initial skepticism about the significance of anti-Glycolipid autoantibodies was replaced by the recognition that these autoantibodies can contribute directly to the pathogenesis of Neuropathy.

GANGLIO Ab PROFILE		
SULFATIDES	→	SENSORY DEMYELINATING POLYNEUROPATHY (IgM)
GQ1b	→	MILLER-FISHER SYNDROME (IgG)
GT1b	→	CANOMAD (IgM)
GT1a	→	CERVICO-BRACHIAL GBS (IgG)
GD3	→	ATAXIC SENSORY NEUROPHATY
GD1b	→	- ACUTE (IgG) - CHRONIC (IgM)
GD1a	→	MOTOR NEUROPATHY (IgM)
GM3	→	CANOMAD (IgM)
GM2	→	CIDP (IgG)
GM1	→	MULTIFOCAL MOTOR NEUROPATHY (IgM) AXONAL FORM OF GBS (AMAN) (IgG)
Funct Control		

2. PRINCIPLE OF THE ASSAY

The ALPHADIA Ganglioprofile kit provides a qualitative in-vitro assay for human autoantibodies (IgG, IgM or Ig G+M) against ten specific Gangliosides: sulfatides, GQ1b, GT1b, GT1a, GD3, GD1b, GD1a, GM3, GM2 and GM1.

During the first incubation time, the specific auto-antibodies present in diluted serum react with the dot ganglioside membrane strips. The strips are then washed to remove unbound antibodies and other serum compounds. In a second incubation time, an anti-human antibody coupled with HRP is added to the tray channel and reacts with the human IgG and IgM antibody immobilized on the membrane strips. After a new washing step to remove the conjugate unbound, specific antibodies are traced by incubation with substrate solution.

3. REAGENTS PROVIDED WITH THE KIT

1. **Coated strips** : 12 or 24 membrane strips coated with ganglioside antigens. Keep unused membrane strips at 4°C, protected from moisture with a desiccant bag.
2. **Positive control** : 1 microtube (30µl) of positive serum
Preservative : NaN₃ (<0.1%)
3. **Enzyme conjugate** : 1 or 2 vials (10 ml) of:
 - anti-hIgG conjugated to horseradish peroxidase (HRPO), for AD-GANG12
 - anti-hIgM conjugated to horseradish peroxidase (HRPO), for AD-GANG12
 - anti-hIg (G+M) conjugated to horseradish peroxidase (HRPO), for AD-GANGS12 and AD-GANGS24Preservative : Neomycin. Ready for use.
4. **Sample diluent** : 1 or 2 vials (32,5 ml) of buffer. Preservative: NaN₃ (<0.1%). Ready for use.
5. **Wash buffer** : 1 vial (40 ml) of phosphate buffer with detergent. Preservative: Thimerosal (<0.1%). Take the vial content to 400 ml (final volume) with distilled water. The diluted washing solution is stable for 1 month at 4°C.
6. **Substrate solution** : 1 or 2 vials (10ml) of tetramethylbenzidine. Ready for use.
7. **Incubation tray** : 1 or 2 x 12 channels

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Adjustable, automatic micropipettes with disposable tips.
- Graduated cylinder.
- Aspiration pump.
- Distilled water.
- Rocking shaker.

5. WARNINGS AND PRECAUTIONS

In order to avoid personal and environmental contamination, the following precautions must be observed:

- Use disposable gloves while handling potentially infectious material and performing the assay.
- Do not pipette reagents by mouth.
- Do not smoke, eat, drink or apply cosmetics during the assay.
- Chromogen should be handled with care. Avoid contact with skin, eyes and mucous membranes. In case of accident rinse thoroughly with running water.
- All material of human origin used for the preparation of this kit tested negative for HBsAg, anti-HIV and anti-HCV. Since no test at present can guarantee complete absence of these viruses, all samples and reagents used for the assay must be considered potentially infectious; therefore, the assay waste must be decontaminated and disposed of, in accordance with established safety procedures. Disposable ignitable material must be incinerated; disposable non-ignitable material must be sterilized in autoclave for at least 1 hour at 121°C.
- Liquid wastes must be added with sodium hypochlorite at a final concentration of 3%. Let the hypochlorite act for at least 30 minutes. Liquid wastes containing acid must be neutralized with appropriate amounts of base before treating with sodium hypochlorite.
- Avoid splashing and aerosol formation; in case of spilling, wash carefully with a 3% sodium hypochlorite solution and dispose of this cleaning liquid as potentially infectious waste.
- Some reagents contain sodium azide as preservative; to prevent

build-up of explosive metal azides in lead and copper plumbing, reagents should be discarded by flushing the drain with large amounts of water.

In order to obtain reproducible results, the following rules must be observed:

- Do not mix reagents from different lot numbers or from other manufacturers
- Strict adherence to the specific time and temperature of incubations is recommended for accurate results
- **Allow to equilibrate the membrane strips to room temperature before open the stored tube**
- **Do not allow the reagents to equilibrate to room temperature before use, keep them at 4°C.**
- Do not use reagents beyond their expiry date
- Incomplete or inefficient washing will cause poor precision and high background
- Use thoroughly clean glassware, free from contamination of metal ions or oxidating substances.
- Use distilled water, stored in clean containers.
- Avoid any contamination among samples; for this purpose, disposable tips should be used for each sample and reagent.
- Microbially contaminated serum or specimens containing heavy, visible particulate should not be used
- Cross contaminations of reagents or sample could cause false results. Use a clean, fresh, disposable pipette tip for each reagent or specimen manipulation
- Do not expose the substrate to light storage or incubation
- Follow exact incubation times.
- A variety of factors influence the assay performance. These include the accuracy and reproducibility of pipetting technique, timing bias during the assay
- **Only dried membrane strip should be evaluated**

6. SPECIMEN COLLECTION

It is recommended to use serum. Highly lipemic or hemolyzed samples must be discarded. Keep samples at 2-8°C for 1 day; for longer periods it is advisable to freeze samples in aliquots at -20°C. Repeated freezing and thawing of samples should be avoided.

7. ASSAY PROCEDURE

ATTENTION :

INCUBATIONS WITH SERUM AND CONJUGATE MUST BE PERFORMED AT 4°C.

To prevent condensation on the strips, open the tubes containing the membrane when the strips have reached the room temperature.

Sample(s) and control predilution : 1:101

- Allocate predilution tubes to sample(s) and control
- Add 500 µl of sample diluent
- Pipette 5 µl of sample(s) and control
- Mix the tubes

Procedure

- Pipette 1 ml of sample diluent, fill each channel containing strips and incubate 10 min **at 4°C** on a rocking shaker
- Pipette 250 µl of prediluted sample and control and add them to the specific channel
- Incubate for 120 minutes **at 4°C** on a rocking shaker
- Eliminate or aspirate the liquid form from each channel
- Wash the membrane strips 3 times for 5min. **at 4°C** with 1ml of diluted washing solution
- Eliminate or aspirate the liquid form from each channel

- Add 750µl of Specific Enzyme Conjugate into the incubation tray.
- Incubate for 60 minutes **at 4°C** on a rocking shaker
- Eliminate or aspirate the liquid form from each channel
- Wash the membrane strips 3 times for 5min. with 1ml of diluted washing solution
- Wash the membrane strips 1 time with 1ml of deionised water
- Eliminate or aspirate the liquid form from each channel
- Add 750µl of substrate-chromogen solution into incubation tray
- Incubate for 10 min. at RT on a rocking shaker
- Eliminate or aspirate the liquid form from each channel
- Stop the reaction with 2ml of deionised water
- Incubate for 5 min. at RT on a rocking shaker
- Remove the strips from the channels.
- Let them air dry for 30 min and read.
- Fix the test strips to the evaluation protocol.

8. ASSAY SCHEME

Incubation tray channel	Control	Sample(s)
Coated strip	1	1
Sample diluent	1000 µl	1000 µl
- Incubate : 10 min at 4°C on a rocking shaker		
Prediluted sample(s)	----	250 µl
Prediluted positive control	250 µl	----
- Incubate : 120 min at 4°C on a rocking shaker		
- Aspirate and wash : 3 x 1000 µl: 5 min at 4°C		
HRPO conjugate	750 µl	750 µl
- Incubate : 60 min at 4°C on a rocking shaker		
- Aspirate and wash : 3 x 1000 µl of washing solution : 5 min at 4°C		
- Aspirate and wash : 1 X 1000 µl of deionised water : 5 min at 4°C		
Substrate solution	750 µl	750 µl
- Incubate : 10 min at RT on a rocking shaker		
- Aspirate the solution		
Desionised water	2000 µl	2000 µl
- Incubate : 5 min at RT on a rocking shaker		
- Air dry		
- Read after 30 min.		

9. INTERPRETATION OF RESULTS

ALPHADOTS strip interpretation

The interpretation of the results consists in a visual examination:

- Any absence of membrane colouring (whereas functional control is reacting) must be considered as a negative result
- Appearance on the membrane of one or several coloured spots (blue deposit with various intensity responses) must be considered as a positive result (+/++/+++).

The intensity of the dot reaction is directly proportional to the amount of specific autoantibodies present in the serum. As antibody levels increase, the dot intensity will also increase, however visual examination alone cannot discern differences between dot intensity that relates as a linear function.

The intensity of the reaction can be appreciated as follows:

Positive	+ /++/+++	A blue coloured distinct dot is easily seen; the result should be interpreted as positive
Weakly positive	±	The dot is not easily seen; the result should be interpreted as borderline
Negative	-	No dot is seen; the result should be interpreted as negative

Weakly positive reaction (±)

With most immunoserology assays, borderline or low intensity reactions indicate minimal antibody levels, only the very high levels in antibodies are significant.

At the final dilution of use (1/510), some spots showing a weak reaction (±) could be difficult to interpret .

This kind of reaction will have thus to be considered as borderline: in this case we advise to re-evaluate the sample at lower dilution (for example final dilution: 1/250) in order to confirm or not its positivity.

However as for any dot blot technique, slight variations of colouring intensities depending on several contributing factors (incubations, T°...) can be observed between different runs : a weakly positive reaction (±) observed at first , may be observed later as negative (-) in another run assay.

Weak reactions may either indicate low-level autoantibodies or non-specific cross-reacting antibodies, which are found in the normal patient population (= natural auto antibodies⁸). In either case, weakly reactive auto antibody (i.e. lower titer) may be reported and should be interpreted with caution since a significant number of normal specimens will have similar reactions.

Quality control

The positive control must be included in each run.

Failure to meet the quality control requirements invalidates the assay.

Positive control: example of result

GANGLIO PROFILE	Positive Control	Interpretation
Sulfatides	●	+
GQ1b		-
GT1b	●	++
GT1a	●	++
GD3	●	+++
GD1b	●	+
GD1a	●	+++
GM3	●	+++
GM2		-
GM1	●	±
Funct Ctl	●	+++

Positive control

- For the interpretation of the result, please refer to the Quality Control sheet included in the kit
- The response observed for positive control can differ from batch to batch
- Positive control reactions may be either slightly darker or lighter depending on incubations conditions.
- **Positive control provided with the kit cannot, in any case, regarded as a cut-off evaluation of the samples**

Functional control

- Functional control present on the bottom of each membrane must show an intense colouring (the intensity of the functional control dot must not be used as a calibrator)
- Functional control assures that the reagents available in the kit are active and that the assay has been performed properly.
- The test should be considered invalid and repeated if the functional control dot is not visible on the strips.

10. CROSS REACTIONS ⁷

	Structure	Also present on:
GM1 – Epitope	A terminal free galactose linked β 1- 3 to N-acetyl - galactosamine	GD1b
GM2 - Epitope	A terminal free N-acetyl – galactosamine linked β 1-4 to galactose	
GM3 - Epitope	A terminal galactose linked β 1-3 to glucose and linked α 2-3 to a sialic acid residue	GD1a, GT1b
GD1a – Epitopes	A central galactose linked α 2-3 to a sialic residue A terminal galactose linked α 2-3 to a sialic acid residue	GM1, GM2, GT1a GT1b
GD1b – Epitopes	A central galactose linked α 2-3 to 2 sialic acid residues linked α 2,8	GT1b, GQ1b
GT1a – Epitope	2 sialic acid residues linked to the terminal galactose	
GQ1b - Epitopes	2 sialic acid residues linked to the terminal galactose 2 sialic acid residues linked to the internal galactose	GT1a GD1b

11. TROUBLE SHOOTING

Unexpected colour development :

- Some sample might give a dark background staining of the membrane. These unspecific reactions have to interpreted as negative
- The wet strips show a dark background colour. This colour vanishes after the strips have dried completely. Only dried strips (30 min) should be evaluated
- Inadequate incubation time and temperature
- Uncontrolled water ingredients
- Substrate or vial used to prepare substrate is contaminated with oxidative active substances

Poor precision :

- non-homogeneous sample after freezing
- turbidity, particles or high lipid content of the sample
- unequal volumes added to the channel
- inadequate elimination of fluids
- washing was incomplete

Don't leave the strips in the channels for drying. Put them on an absorbent tissue.

12. LIMITATIONS OF THE PROCEDURE

- All ALPHADIA kits are a diagnostic aid and by themselves are not diagnostic. The data must be interpreted in conjunction with the clinical evaluation or/and other laboratory findings.
- The kit should not be used beyond the expiration date on the kit label.
- The procedure must be carry out at 4°C. A modification of this temperature induce a loss of sensibility of the dosage.

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1. Fredman, P. and A.Lekman. 1997. Glycosphingolipids as potential diagnostic markers and/or antigens in neurological disorders. *Neurochem. Res.* 22:1071-1083.
2. Willison, H.J. and N.Yuki. 2002. Peripheral neuropathies and anti-glycolipid antibodies. *Brain* 125:2591-2625.
3. Vincent, A., M.Roberts, H.Willison, B.Lang, and J.Newsom-Davis. 1995. Autoantibodies, neurotoxins and the nervous system. *J. Physiol Paris* 89:129-136.
4. Fredman, P. 1998. The role of antiglycolipid antibodies in neurological disorders. *Ann. N. Y. Acad. Sci.* 845:341-352.
5. Goffette, S., A.Jeanjean, F.Pierret, A.Peeters, and C.J.Sindic. 1998. Clinical relevance of the determination of anti-GQ1b antibodies in Miller Fisher and Guillain-Barre syndromes. *Acta Neurol. Belg.* 98:322-326.
6. Alaedini, A., H.W.Sander, A.P.Hays, and N.Latov. 2003. Antiganglioside antibodies in multifocal acquired sensory and motor neuropathy. *Arch. Neurol.* 60:42-46.
7. Humbel R.L., C.H.L. Luxembourg.
8. Caudie C, Vial C, Bancel J, Petiot P, Antoine JC, Gonnaud PM; Measurement of antiganglioside autoantibodies by immunodot-blot assay: clinical importance in peripheral neuropathies. *Ann Biol Clin.* 1999 Sep-Oct;57(5):579-88.