ENZYWELL

HELICOBACTER PYLORI IgG

REF 91060 (96 tests)

Manufactured by: DIESSE Diagnostica Senese
Via delle Rose 10
53035 Monteriggioni (Siena) - Italy

INDICE / INDEX

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(english)

1. INTENDED USE:

IMMUNOENZYMATIC KIT FOR THE QUALITATIVE AND QUANTITATIVE DETERMINATION OF IgG-CLASS ANTIBODIES TO HELICOBACTER PYLORI IN HUMAN SERUM

2. SUMMARY AND EXPLANATION OF THE TEST:

In 1983, Warren and Marshall identified *Helicobacter pylori*, a new gram-negative bacterial pathogen, in patients suffering from gastritis, and this finding led to studies on the relationship between bacterial infection and chronic gastric disease. The pathogen has been shown to be associated with peptic ulcer, chronic gastritis type B and duodenitis. It has been demonstrated that in patients with gastritis, eradication of the bacteria led to healing of the anatomical lesion.

Diagnostic procedures for the detection of the organism generally involve invasive (gastroscopic) techniques for sample collection.

However, a specific immune response is seen in infected patients. The serological test thus represents a useful alternative to the invasive bioptic technique. IgG levels rise with infection and remain constantly high until the infection is eliminated. The efficacy of antimicrobial therapy can therefore be monitored via changes in specific IgG antibody.

3. PRINCIPLE OF THE TEST:

The test is based on the ELISA technique (Enzyme linked Immunosorbent Assay).

The antigen is bound to the solid phase. The specific immunoglobulins are bound to the antigen through incubation with dilute human serum.

After washings to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of human IgG monoclonal antibodies labelled with peroxidase.

The unbound conjugate is eliminated, and the peroxidase substrate added.

The density of the colour which develops is proportional to the concentration of specific antibodies present in the serum sample. When the enzymatic reaction is interrupted by the addition of a sulfuric acid solution, the resulting yellow colour can be easily read on an ELISA microplate reader.

4. KIT CONTENTS AND REAGENT PREPARATION:

- Reagents are sufficient for 96 determinations.

- **Bring reagents to room temperature before use.**

MT PLATE MICROPLATE. 12x8 wells coated with *Helicobacter pylori* antigens.

Use: open the package at the opposite end from the code (P, followed by the lot number) which serves for identification purposes, remove the support and strips to be used from the foil package, and place the unused strips in the polythene bag with the silica gel, expel the air and seal by pressing the closure.

CAL CALIBRATORS. 5 x 1.6 mL.

Contents: Diluted human serum at known concentrations of antibodies to *Helicobacter pylori*, in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%. Liquid, ready for use without further dilution.

The 5 calibrators have the following titers expressed in Arbitrary Units/ml:

Cal 1 5 AU/ml (Negative Calibrator)

Cal 2 (Pediatric Cut-off) 10 AU/ml

Cal 3 (Adult Cut-off) 15 AU/ml

Cal 4 50 AU/ml

Cal 5 100 AU/ml

Colour: the colour of the calibrators is proportional to the relative antibody titer.

CONJ CONJUGATE. 1x16 mL.

Contents: monoclonal antibodies anti-human IgG labelled with Peroxidase, in phosphate buffer with phenol 0.05% and Bronidox 0.02%. Ready for use without further dilution.

WASH BUF WASH BUFFER 10X (PF93603). 1 x 100 mL. **INTERCHANGEABLE BETWEEN LOTS**

Contents: Phosphate buffered saline, concentrated 10 times; contains Brij 0.5% .

Preparation: dilute the required volume 1:10 with distilled water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

SAMP DIL 50x DILUENT 50X (PF93601). 1x4.5 mL. For dilution of serum samples. **INTERCHANGEABLE BETWEEN LOTS**

Contents: Proteic solution concentrated 50 times, with added phenol 0.05% and Bronidox 0.02%.

Preparation: Dilute the required volume 1:50 in the washing buffer to obtain the diluent ready for use.

SUBS TMB SUBSTRATE (PF93619). 12 mL. Ready for use. **INTERCHANGEABLE BETWEEN LOTS**

Contents: Tetramethylbenzidine 0.26 mg/mL and hydrogen peroxide 0.01% stabilised in citrate buffer 0.05 mol/L (pH 3.8).

H₂SO₄ 0.3 M STOP SOLUTION (PF93602). 1x16 mL. **INTERCHANGEABLE BETWEEN LOTS**
H₂SO₄ 0.3 mol/L, in solution ready for use.

ADHESIVE FILMS (2).

POLYTHENE BAG (1).

MATERIALS REQUIRED BUT NOT PROVIDED.

- Incubator at 37°C
- Microplate reader (wave length 450 or 450/620 nm, with linearity up to OD >= 2000)
- Microplate washer (preferable) able to dispense volumes ranging between 225-375 µl
- Distilled or deionised water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µl solution
- Disposable gloves
- Timer
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials
- Absorbent tissue.

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation

REAGENT	CONDITIONS
Microplate	6 weeks at 2/8°C, polythene bag
Calibrators	6 weeks at 2/8°C
Conjugate	6 weeks at 2/8°C
Substrate	up to the expiry date at 2/8°C, 1 week at 15-30°C; store in the dark
Sample Diluent	ready for use, 2 weeks at 2/8°C
Wash Buffer	2 weeks at 2/8°C, 5 days at 15/30°C.
Stop Solution	up to the expiry date at 2/8°C

6. PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY.

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The Wash Buffer contains detergents
 - b) The conjugate contains phenol
 - c) The substrate is acid
 - d) The controls contain 0.09% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides; dilute with large amounts of water to eliminate.
If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Non-disposable apparatus should be sterilized after use. The preferred method is to autoclave for 1 h at 121°C; disposables should be autoclaved or incinerated.
4. Sulphuric acid required for the Stop Solution and hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If they come into contact with the skin or eyes, wash thoroughly with water.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature. **It is important to work at the correct temperature. Check that the thermostat does not go below 35°C or over 39°C.** Open the envelope containing the strips after at least ½ hr at room temperature.
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
4. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionized water.
5. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
6. Do not allow wells to become dry during the assay procedure.
7. Care must be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for exclusive use with the various reagents.
8. Care should be taken to avoid touching or splashing the rim of the well with conjugate. Do not "blow-out" from microplates.
9. Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimized by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at 37°C either in a water bath with a rack or float to support the plates if necessary, or in an incubator. Alternatively, plates can be incubated in an approved analyzer. See the appropriate operating manual for further details. CO₂ incubators must not be used.
10. Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate.
11. Use of highly hemolyzed samples, incompletely clotted sera, or samples with microbial contamination may give rise to erroneous results.
12. Use of the kit with automatic instruments must be validated by the user.
13. For each instrument used, read the manufacturer's instructions manual carefully to obtain additional information on the following points:
 - installation and particular requisites
 - operating principles, instructions, precautions and risks
 - manufacturer's specifications and instrument performance
- servicing and maintenance.
-

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at

-20°C, and can be thawed a maximum of 3 times. Defrosted samples must be carefully mixed before performing the test. Heat inactivation can lead to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric or contaminated samples should be avoided. If a new sample cannot be obtained, such samples should be clarified by filtration (0.45 µm) or centrifugation (3000 rpm x 10').

The test is not applicable to human plasma.

8. TEST PROCEDURE:

- Prepare the required number of strips.
- Prepare the washing buffer by diluting the Wash Buffer 10x (100 mL + 900 mL H₂O).
- Prepare the required diluent for samples by adding 1 part of Diluent 50x to 49 parts of the diluted washing buffer (example: 2 mL + 98 mL of diluted washing buffer).

Dilute the samples 1:101 distributing 10 µL of serum into 1 mL of diluent. Leave one well empty for the blank (only the substrate will be used). Pipet 100 µL of UNDILUTED calibrators in the wells of a strip, preferably in duplicate, and 100 µL of diluted sample in the remaining wells.

Cover the wells with protective film and incubate for 45 minutes at 37°C. After washing four times (300 µL), soaking the wells in the washing solution for 30 seconds for each cycle, add 100 µL of conjugate to each well and incubate again for 45 minutes at 37°C, covering the wells with the protective film. Wash the plate again 4 times, as described above. Finally, add the Substrate (100 µL/well).

After 15 minutes at room temperature, stop the enzymatic reaction with 100 µL of Stop Solution.

Read the optical densities (O.D.) at 450 nm or at 450/620 nm within 30 min. Take a new reading at 405 nm if the O.D. are higher than 2000.

9. TEST PROCEDURE FOR HELICOBACTER PYLORI IgG

- | | |
|--------|--|
| STEP 1 | Place 100 µL of diluted sample//controls in the wells of the microplate.
<input type="checkbox"/> |
| | Incubate for 45 min. at 37°C
<input type="checkbox"/> |
| | Wash 4 times (300 µL)
<input type="checkbox"/> |
| STEP 2 | Add 100 µL of conjugate to each well
<input type="checkbox"/> |
| | Incubate for 45 min. at 37°C
<input type="checkbox"/> |
| | Wash 4 times (300 µL)
<input type="checkbox"/> |
| STEP 3 | Add 100 µL of Substrate to each well
<input type="checkbox"/> |
| | Incubate for 15 min. at R.T.
<input type="checkbox"/> |
| STEP 4 | Add 100 µL of Stop Solution
<input type="checkbox"/> |
| | Read absorbance at 450 nm within 30 min |

10. TEST VALIDATION

Subtract the value of the optical density (O.D.) of the blank (≤ 0.150) from all the other readings. The O.D. of Calibrator 3 must be ≥ 0.2 . Calibrator 5 must have an optical density (O.D.) at least 1.5 times that of Calibrator 3.

11. INTERPRETATION OF THE RESULTS

Qualitative results

The Cut-off for adult sera corresponds to Calibrator 3. The Cut-off for pediatric sera corresponds to Calibrator 2.

Calculate the ratio between the average O.D. value of the sample and that of the Cut-off (INDEX). The sample is considered:

Positive: if the ratio is > 1.2 .

Doubtful: $\pm 20\%$ of the Cut-Off.

Negative: if the ratio is < 0.8 .

If the result is doubtful, repeat the test. If it remains doubtful, collect a new serum sample.

Quantitative results

Report the OD of the calibrators on a graph after subtracting the OD of the blank. The corresponding titer of the test sample can be obtained by extrapolation.

NOTE: A standard curve must be performed for each run.

If the O.D. of any sample or calibrator is over 2.0, take the reading at 405 nm and multiply the value by 3.

The degree of immunity can be interpreted as follows:

IMMUNE: when the anti-H. Pylori IgG concentration in the sample is > 10 AU/ml (for pediatric patients)
 > 15 AU/ml (for adult patients)

NON IMMUNE: when the anti-H. pylori IgG concentration is < 5 AU/mL

DOUBTFUL : if the result is between the two values. In this case it is advisable to repeat the test in duplicate.

12. LIMITATIONS OF THE PROCEDURE

A serum sample obtained during the early acute phase of infection may be negative by this procedure.

Highly lipemic, icteric, hemolysed or heat-inactivated sera may cause false results and should not be used.

As with all serological tests, the results should be interpreted in conjunction with other clinical and diagnostic data.

13. ANALYTICAL SPECIFICITY

The analytical specificity of the Enzywell H. Pylori IgG kit was evaluated in an inhibition test. Campylobacter coli, Campylobacter fetus, Campylobacter jejuni and E. coli were added to patients' serum samples. No differences were observed between the reference samples and those containing the antigens reported above. It can therefore be concluded that the test is specific for anti-Helicobacter pylori IgG.

14. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

In a clinical trial carried out in a hospital laboratory, 199 human sera were analyzed, 58 of which came from the Gastroenterology Department. Out of the total, 138 gave a positive result.

The samples were analyzed in Western blot, considered the reference method.

Results

		REFERENCE	
		+	-
DIESSE	+	132	6
	-	1	60

The Helicobacter pylori kit result 6 false negative and 1 false positive.

The sensitivity and specificity were respectively 95.7% and 98.4%.

15. PRECISION

Within run Precision on three different lots

Cut Off $n=16$	Batch. N. 78	Batch. N. 79	Batch. N. 80
O.D.	0.419	0.608	0.6
CV%	12	15	10

Between run Precision

Sample	AU/ml				Average	CV%
	I run	II run	III run	Average		
PYG1	6.2	7	6.1	6.4	8	

PYG2	24.8	26.1	23.9	24.9	4
PYG3	81.3	88.4	84.8	84.8	4

16. TROUBLE SHOOTING GUIDE

PROBLEM	POSSIBLE SOURCE	TEST OR ACTION
Invalid run (all negative)	One or more reagents not added or added in wrong sequence	Recheck procedure Check for unused solutions. Repeat test.
	Unreactive plate	Check the code on the package containing the plate (see instructions for use point 4 for correct code).
		Check for moisture in unused plate. (Silica gel dessiccant must be pale yellow).Repeat test
Invalid run (all positive)	Contamination of substrate	Take new aliquot of substrate.
	Inadequate washing	Ensure that wash apparatus works well
Poor precision	Incomplete washing of wells	Ensure that wash apparatus works well
	Inadequate aspiration of wells	Ensure that wash apparatus works well
	Pipetting error	Check pipette function
	Reagent addition too slow	Avoid drying of the plate after washing step. Add reagents immediately
	Presence of bubbles	Avoid air bubbles during pipetting.
	Optical pathway not clean	Check instrument light source and detector for dirt. Wipe bottom of plate with soft tissue.
Inadequate Color development	Incorrect incubation times or temperature	Check for temperature control and time monitoring
		Adhere to recommended instruction for use.
	Inadequate volume of substrate added to the plate	Check pipette function.

17. BIBLIOGRAFIA / REFERENCES

1. Marshall B.J. and Warren J.R.: Unidentified curved bacilli in the stomach of patients with gastritis and peptic ulceration. *Lancet* i, 1311 (1984).
2. Jones D.M., Lessels A.M., Eldridge J.: Campylobacter-like organisms on the gastric mucosa: culture, histological and serological studies. *J. Clin. Pathol.* 37: 1002 (1984).
3. Blaser M.J.: H. pylori and the pathogenesis of gastroduodenal inflammation. *J. Inf. Dis.* 161: 626 (1990).
4. Valle J., Seppälä, K., Sipponen P., Kasunen T.U.: Disappearance of gastritis after eradication of H. pylori: a morphometric study. *Scand. J. Gastroenterology* 26: 1057 (1991).
5. G.B. Wisdom: Enzyme-Immunoassay. *Clin. Chem.* 22: 1243 (1976).

DIESSE Diagnostica Senese
Via delle Rose 10
53035 Monteriggini (Siena), Italy.