ENZYMATIC HOMOCYSTEINE ASSAY
DOSAGE ENZYMATIQUE POUR L' HOMOCYSTEINE
ENZYMATISCHER HOMOCYSTEINASSAY
ANALISI OMOCISTEINA ENZIMATICA
ENSAYO DE HOMOCISTEÍNA ENZIMÁTICA
ENZYMATISK HOMOCYSTEINANALYS
ENZYMATISK HOMOCYSTEINANALYSE
TEST NA ENZYMATICKÝ HOMOCYSTEIN
ENZYMATICKÝ KVANTITATÍVNY ROZBOR
HOMOCYSTEÍNU

IVD

REF FHER100

For professional use only
Usage réservé aux professionnels
Nur für den Fachgebrauch
Solo per uso professionale
Para uso profesional únicamente
Endast för professionellt bruk
Kun til fagligt brug
Pouze pro odborné použití
Iba na profesionálne použitie

Axis-Shield Diagnostics Limited
The Technology Park, Dundee, DD2 1XA, United Kingdom
Tel: +44 (0) 1382 422000 Fax: +44 (0) 1382 422088
E-mail: shield@uk.axis-shield.com
Web: www.axis-shield.com

ENGLISH:**INTENDED USE**




The Enzymatic Homocysteine Assay is intended for *in vitro* quantitative determination of total homocysteine in human serum and plasma. This assay is for professional use only.

INTRODUCTION

Homocysteine (HCY) is a by-product of protein metabolism, specifically the breakdown of the essential amino acid, methionine. It is not by itself incorporated into proteins but is further broken down by enzymes into cysteine or is converted back to methionine. In small amounts homocysteine is not harmful to the body or blood vessels, but when excess amounts accumulate in the blood stream, arterial vessels may be damaged and the resulting inflammation may eventually cause blockage of blood to the heart. Clinical studies began in 1969, when Dr. Kilmer S. McCully reported that children born with a genetic error of metabolism that causes elevated homocysteine levels (homocystinuria) died with advanced artery disease at a very young age¹. More recent studies are yielding evidence that elevated blood levels of homocysteine have a predictive value for risk of coronary artery disease similar to that of elevated cholesterol levels^{2,3}. Recent evidence has also implicated elevated blood levels of homocysteine in miscarriages and birth defects⁴.

PRINCIPLE OF THE ASSAY

Bound or dimerised homocysteine (oxidised form) is reduced to free homocysteine, which then reacts with serine catalysed by cystathionine beta-synthase (CBS) to form L-cystathionine. Cystathionine in turn is broken down by cystathionine beta-lyase (CBL) to form homocysteine, pyruvate and ammonia. Pyruvate is then converted by lactate dehydrogenase (LDH) to lactate with NADH as coenzyme. The rate of NADH conversion to NAD is directly proportional to the concentration of homocysteine ($\Delta A_{340\text{ nm}}$).

**KIT COMPONENTS**

Enzymatic Homocysteine Reagent R1	Clear odourless liquid	1 x 36ml in Amber Vial	NADH (0.56mM), LDH (65KU/L), Serine (1.3mM), Trizma Base 1-10%, Trizma Hydrochloride 1-10%, Sodium Azide <1%. Ready-to-use	
Enzymatic Homocysteine Reagent R2	Clear odourless liquid	1 x 15ml in Clear Vial	Reductant (TCEP;25mM). Ready-to-use	
Enzymatic Homocysteine Reagent R3	Pale yellow odourless liquid	1 x 5ml in Amber Vial	Cycling Enzymes CBS (0.748KU/L) and CBL (16.4KU/L) Sodium Azide <1%. Ready-to-use	
Enzymatic Homocysteine Calibrator 1	Clear odourless liquid	1 x 3ml in Clear Vial (Blue Cap)	Aqueous homocysteine blank (0µmol/L). Ready-to-use	
Enzymatic Homocysteine Calibrator 2	Clear odourless liquid	1 x 3ml in Clear Vial (Red Cap)	Aqueous homocysteine solution (27µmol/L). Ready-to-use	

The calibrators are prepared gravimetrically and are traceable to a designated measurement procedure (HPLC). The values are printed on the labels (0µmol/L and 27µmol/L).

STORAGE OF REAGENTS

Handling and Procedural Notes

The Enzymatic Homocysteine Reagents and Calibrators should be stored refrigerated at 2-8°C. **DO NOT FREEZE.** The reagents and calibrators are stable when stored as supplied until the expiration dates on the labels. The reagents are stable on board the instrument for up to 30 days. Do not mix lots of reagents.

Specimen Collection and Handling

Fresh EDTA or heparinised plasma free of haemolysis or turbidity is the recommended specimen for measurement of homocysteine. Serum collected in serum separator tubes may also be used. However, it is not recommended to use individual patient results from serum, heparinised plasma and EDTA plasma interchangeably. All samples must be placed on ice immediately after collection and should be centrifuged and separated from cells as soon as possible (within one hour). Severely lipaemic specimens are unsatisfactory. Samples should be stored tightly capped at 2-8°C for up to 48 hours, or frozen at -20°C if testing is delayed⁶.

WARNINGS AND PRECAUTIONS

For in vitro diagnostic use only

Safety Precautions

1. Adhere strictly to the instructions in this leaflet, particularly for handling and storage conditions.
2. Reagents R1 and R3 contain sodium azide which can react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with large quantities of water to prevent azide build-up.
3. Material safety data sheets for all hazardous components contained in this kit are available upon request from Axis-Shield Diagnostics Ltd.

Reagent R1
Reagent R3
Harmful

R22: Harmful if swallowed.

R32: Contact with acids liberates very toxic gas.

S36/37/39: Wear suitable protective clothing, gloves and eye/face protection.

S29/35: Do not empty into drains; dispose of this material and its container in a safe way.

S46: If swallowed, seek medical advice immediately and show this container or label.

PREPARATION

Materials Required but not Provided

1. Deionised Water.
2. An analyser capable of dispensing three reagents and of measuring absorbance at 340nm with temperature control (37°C).

ASSAY PROTOCOL⁸

RESULTS

Results are printed out by the Cobas MIRA in µmol/L.

QUALITY CONTROL

Ensure that adequate maintenance and calibration is performed according to the manufacturer's instructions.

Assayed control materials with values for homocysteine in both the normal and abnormal ranges should be tested to validate reagent performance. The range of acceptable control limits should be established by individual laboratories. Axis-Shield provide a set of low, medium and high controls.

REF FHER200

The controls contain L-homocysteine in processed human serum at the following concentrations:

<i>Control</i>	<i>Mean Value HCY (µmol/L)</i>	<i>Range HCY (µmol/L)</i>
<i>Control L</i>	7.0	5.6 – 8.4
<i>Control M</i>	12.5	10.0 – 15.0
<i>Control H</i>	25.0	20.0 – 30.0

PERFORMANCE DATA

(Established on the Cobas MIRA)

Accuracy

Correlation studies were performed with fresh plasma from adults. All samples were analysed using the enzymatic assay, Abbott FPIA assay and the Pacific Biometrics, Inc. HPLC assay. The data obtained gave the following statistical values:

Comparison Method	Abbott FPIA	HPLC
Number of samples	85	30
Slope of regression line	0.964	1.059
Y-Intercept	0.728	-0.035
Correlation coefficient	0.993	0.989
Range ($\mu\text{mol/L}$)	4.6-51.8	5.1-26.0

Precision

Precision studies were conducted utilising the Enzymatic Homocysteine Assay according to the NCCLS Document EP5-A⁷ protocol. Each human EDTA plasma sample was run in duplicate on forty runs spread over 20 days. Results are summarised below:

Sample	Mean ($\mu\text{mol/L}$)	Within Run CV%	Between Run CV%	Total CV%
1	4.6	6.6	3.1	7.3
2	11.6	3.5	0.7	4.0
3	27.0	1.6	1.5	2.6

Analytical Sensitivity

The analytical sensitivity of the Enzymatic Homocysteine assay was found to be 0.40 $\mu\text{mol/L}$.

Analytical Specificity

The specificity of the assay was assessed for the interfering substances listed in the table below:

Interfering Substance	Interfering Substance Concentration	% Interference
Bilirubin	20 mg/dL	< 5%
Haemoglobin	1000 mg/dL	< 5%
Lipaemia	1000 mg/dL	< 10%
Glutathione	20 $\mu\text{mol/L}$	< 5%
Methionine	800 $\mu\text{mol/L}$	< 5%
Cysteine	200 $\mu\text{mol/L}$	< 5%
Pyruvate	1000 $\mu\text{mol/L}$	< 5%

None of these substances interfered significantly (<10%) in the assay.

Refer to references (5,6) for possible interferences caused by drugs, disease or preanalytical variables

Carryover

Carryover studies show that carryover is less than the sensitivity of the assay.

EXPECTED VALUES

A study was conducted using samples from 131 adults (59 males and 72 females, ages 18-86) who were seemingly healthy. Plasma (EDTA) samples were assayed for total homocysteine by the Enzymatic Homocysteine Assay and the reference range for homocysteine by this method is:

Adult 4.0-15.4 $\mu\text{mol/L}$

These values are to be considered as guidelines for expected values. Each laboratory is encouraged to establish a range of normal values for the population they serve.

LIMITATIONS OF USE

1. The linearity of the Enzymatic Homocysteine Assay when run as directed is 50 $\mu\text{mol/L}$. Samples with values greater than 50 $\mu\text{mol/L}$ should be diluted 1 in 3 in saline and the result multiplied by three ($\times 3$).
2. The R1 reagent should be clear. It should be discarded if it becomes turbid or the initial absorbance is less than 0.700.
3. Cystathionine is measured with Homocysteine, but in the general population the cystathionine level (0.065 to 0.3 $\mu\text{mol/L}$) has a negligible effect. In some occasions the cystathionine levels may rise in renal disease patients with severe metabolic disturbances and cause greater than 20% interference. Although some studies have shown that no patients are affected, others have shown that this effect may be seen in up to 10% of renal disease patients.
4. **Drugs, such as methotrexate, carbamazepine, phenytoin, nitrous oxide, or 6 azauridine triacetate may affect the plasma homocysteine concentration^{5,6}.**

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