Afinion™ HbA1c

REF 1115016

Hemoglobin A1c test	GB	5...12
Hämoglobin-A1c-Test	DE	13...21
Dosage de l'hémoglobine A1c	FR	22...30
Test per l'emoglobina A1c	IT	31...39
Test HbA1c	ES	40...48
Teste Hemoglobina A1c	PT	49...57
Hemoglobine A1c test	NL	58...66
Αιμοσφαιρίνη A1c δοκιμή	GR	67...76

CE

afinion™
POINT OF CARE CENTER

SYMBOLS/ SYMBOLE/ SYMBOLES/ SIMBOLI / SÍMBOLOS/ SÍMBOLOS/ SYMBOLLEN/ ΣΥΜΒΟΛΑ

CE

Conformity to the European directive 98/79/EC on *in vitro* diagnostic medical devices
 Konformität mit der europäischen Richtlinie 98/79/EC über medizinische Geräte zur *In-vitro*-Diagnostik
 Conformité avec la directive européenne 98/79/CE sur les appareils médicaux de diagnostic *in vitro*
 Conformità con la direttiva europea 98/79/CE sui dispositivi medicali per la diagnostica *in vitro*
 Conformidad con la directiva europea 98/79/EC sobre dispositivos médicos para diagnóstico *in vitro*
 Conformidade com a Directiva Europeia 98/79/EC de dispositivos médicos de diagnóstico *in vitro*
 Voldoet aan de Europese richtlijn 98/79/EG betreffende medische hulpmiddelen voor *in-vitro*diagnostiek
 Συμμορφώνεται με την ευρωπαϊκή οδηγία 98/79/ΕΚ περί *in vitro* διαγνωστικών ιατρικών συσκευών

IVD

In Vitro Diagnostic Medical Device
 Medizinisches Gerät zur *In-vitro*-Diagnostik
 Dispositif médical de diagnostic *in vitro*
 Dispositivo medico per la diagnostica *in vitro*
 Dispositivo médico para diagnóstico *in vitro*
 Dispositivo médico de diagnóstico *in vitro*
In-vitro-diagnostisch medisch apparaat
In Vitro Διαγνωστική Ιατρική Συσκευή

REF

Catalogue number
 Bestellnummer
 Numéro de référence du catalogue
 Numero di catalogo
 Número de catálogo
 Número do catálogo
 Catalogusnummer
 Αριθμός καταλόγου

LOT

Lot number
 Chargen-Nr.
 Numéro de lot
 Numero di lotto
 Número de lote
 Número de lote
 Partijnummer
 Αριθμός Παρτίδας

TEST CARTRIDGE

Test Cartridge
 Testkassette
 Cassette de dosage
 Cartuccia per test
 Cartucho de análisis
 Cartucho de ensaio
 Testcartridge
 Κασέτες Εξέτασης

Contents sufficient for "Σ" number of tests
 Inhalt ausreichend für "Σ" Tests
 Contenu suffisant pour "Σ" dosages
 Contenuto sufficiente per un numero "Σ" di test
 Contenido suficiente para "Σ" cantidad de pruebas
 Contém suficiente para "Σ" testes
 Inhoud voldoende voor "Σ" aantal testen
 Επαρκή περιεχόμενα για "Σ" αριθμό εξετάσεων


Consult the Afinion™ user instructions
 Bedienungsanleitungen für Afinion™ beachten
 Consultez les instructions d'utilisation Afinion™
 Consultare le istruzioni per l' utilizzo dell' Afinion™
 Consulte las instrucciones de uso de Afinion™
 Consulte as instruções para o utilizador da Afinion™
 Zie de gebruikersinstructies voor de Afinion™
 Συμβουλευτείτε τις οδηγίες χρήσης του Afinion™

Caution, consult instructions for use
 Achtung, Bedienungsanleitung beachten
 Attention, consultez les instructions d'utilisation
 Attenzione, consultare le istruzioni d' uso
 Precaución, consulte las instrucciones de uso
 Atenção: consulte as instruções de utilização
 Opgelet, raadpleeg de gebruiksaanwijzing
 Προσοχή, συμβουλευτείτε τις οδηγίες χρήσης

Expiry date (year-month)
 Verfallsdatum (Jahr-Monat)
 Date de péremption (année-mois)
 Data di scadenza (anno-mese)
 Fecha de caducidad (año-mes)
 Prazo de validade (ano-mês)
 Uiterste gebruiksdatum (jaar – maand)
 Ημερομηνία λήξης (έτος-μήνας)

Storage temperature 2-8°C (36-46°F)
 Lagertemperatur 2-8°C (36-46°F)
 Température de conservation 2-8°C (36-46°F)
 Temperatura di conservazione 2-8°C (36-46°F)
 Temperatura de almacenamiento 2-8°C (36-46°F)
 Temperatura de conservação 2°C-8°C (36-46°F)
 Bewaartemperatuur 2-8°C (36-46°F)
 Θερμοκρασία φύλαξης 2-8°C (36-46°F)

Manufacturer
 Hersteller
 Fabricant
 Produttore
 Fabricante
 Fabricante
 Fabrikant
 Κατασκευαστής

For use with the Afinion™ AS100 Analyzer.

PRODUCT DESCRIPTION

Intended use

Afinion™ HbA1c is an *in vitro* diagnostic test for quantitative determination of glycated hemoglobin (hemoglobin A1c, HbA1c) in human whole blood. The measure of HbA1c is recommended as a marker of long-term metabolic control in persons with diabetes mellitus¹.

Summary and explanation of the test

The human erythrocyte is freely permeable to glucose. Within each erythrocyte a slow, continuous, non-enzymatic process between hemoglobin A and various sugars takes place. The product formed is known as glycated hemoglobin, or glycohemoglobin².

The chronic elevated blood sugar level of persons with diabetes mellitus will over time cause damage to the small vessels of the body. This damage develops slowly over years and is known to cause late complications³.

Good metabolic control, i.e. lowering the HbA1c concentration, has proven to delay the onset and slow the progression of diabetes late complications^{3,4}.

Principle of the assay

Afinion™ HbA1c is a fully automated boronate affinity assay for the determination of hemoglobin A1c in human whole blood.

The Afinion™ HbA1c Test Cartridge contains all the reagents necessary for the determination of the HbA1c concentration. The sample material is collected using the integrated sampling device and the Test Cartridge is placed in the Afinion™ AS100 Analyzer. The blood sample is then automatically diluted and mixed with a liquid that releases hemoglobin from the erythrocytes. The hemoglobin precipitates. This sample mixture is transferred to a blue boronic acid conjugate, which binds to the cis-diols of glycated hemoglobin. This reaction mixture is soaked through a filter membrane and all precipitated hemoglobin, conjugate-bound and unbound (i.e. glycated and non-glycated hemoglobin) remains on the membrane. Any excess of conjugate is removed with a washing reagent.

The Afinion™ AS100 Analyzer evaluates the precipitate on the membrane. By measuring the reflectance, the blue (glycated hemoglobin) and the red (total hemoglobin) colour intensities are evaluated, the ratio between them being proportional to the percentage of HbA1c in the sample. The HbA1c concentration is displayed on the Afinion™ AS100 Analyzer in the units mmol/mol, percentage (%), as the calculated estimated average glucose (eAG) or in combinations of these.

Kit contents (per 15 tests unit)


- 15 Test Cartridges packaged separately in foil pouches with a desiccant bag
- 1 Package Insert

Materials required but not provided with the kit

- Afinion™ AS100 Analyzer
- Afinion™ HbA1c Control
- Standard blood collection equipment

Description of the Test Cartridge

The main components of the Test Cartridge are the sampling device and the reagent container. The Test Cartridge has a handle, a barcode label with lot-specific information and an area for sample ID. See figure and table below.

Component	Function/composition
1 Sampling device a. Closed position b. Lifted position	For collection of patient sample or control.
2 Capillary	1.5 µL glass capillary to be filled with sample material.
3 Reaction wells a. Conjugate b. Membrane tube c. Washing solution d. Reconstitution reagent e. Empty	Contains reagents necessary for one test: Patented blue boronic acid conjugate. Tube with a polyethersulfone membrane. Morpholine buffered sodium chloride with detergents and preservative. HEPES buffered sodium chloride with lysis and precipitation agents. N/A
4 Handle	For correct finger grip.
5 Barcode label	Contains assay- and lot-specific information for the Analyzer.
6 Optical reading area	Area for transmission measurement.
7 ID area	Space for written or labelled sample identification.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use.
- Do not use Test Cartridges after the expiry date or if the Test Cartridges have not been stored in accordance with recommendations.
- Do not use the Test Cartridge if the foil pouch or the Test Cartridge has been damaged.

- Each foil pouch contains a desiccant bag with 1 g silica gel. This material shall not be used in the assay. Discard the desiccant bag in a suitable container. Do not swallow.
- Do not use the Test Cartridge if the desiccant bag is damaged and desiccant particles are found on the Test Cartridge.
- The Afinion™ HbA1c Test Cartridge contains sodium azide as a preservative. The concentration is < 0.1%, which is below that which is considered hazardous in normal use⁶. In case of leakage, avoid contact with eyes and skin. Wash with plenty of water.
- Do not re-use any part of the Test Cartridge.
- The used Test Cartridges, sampling equipment, patient samples and controls are potentially infectious. The Test Cartridges should be disposed of immediately after use. Proper handling and disposal methods should be followed in accordance with local or national regulations. Use gloves.

STORAGE

Refrigerated storage (2-8°C)

- The Afinion™ HbA1c Test Cartridges are stable until the expiry date only when stored refrigerated in sealed foil pouches. The expiry date is the last day of the month stated on the foil pouch and kit container.
- Do not freeze.

Room temperature storage (15-25°C)

- The Afinion™ HbA1c Test Cartridges can be stored in unopened foil pouches at room temperature for 90 days. Note the date of removal from the refrigerator on the kit container.
- Avoid exposure to direct sunlight.

Opened foil pouch

- The Test Cartridge must be used within 10 minutes after opening the foil pouch.
- Avoid exposure to direct sunlight.
- Avoid relative humidity above 90% (non-condensing).

SAMPLE MATERIAL

The following sample materials can be used with the Afinion™ HbA1c test:

- Capillary blood (from finger prick)
- Venous whole blood with anticoagulants (EDTA, heparin, citrate or NaF)
- Afinion™ HbA1c Control

Sample storage

- Capillary blood without anticoagulants cannot be stored.
- Venous whole blood with anticoagulants (EDTA, heparin or citrate) may be stored refrigerated for 10 days. Do not freeze.
- Consult the Afinion™ HbA1c Control Package Insert for storage of control materials.

PREPARING FOR ANALYSIS

 Consult the Afinion™ AS100 Analyzer User Manual for detailed instructions on how to analyse a patient or control sample. The Afinion™ HbA1c Quick Guide also provides an illustrated step-by-step procedure.

- The Afinion™ HbA1c Test Cartridge must reach an operating temperature of 18-30°C before use. Upon removal from refrigerated storage, leave the Test Cartridge in unopened foil pouch for at least 15 minutes.
- Open the foil pouch just before use.
- Do not touch the Test Cartridge optical reading area.
- Label the Test Cartridge with patient or control ID. Use the dedicated ID area.

COLLECTING A SAMPLE

Capillary blood

- Make sure that the finger is clean, warm and dry.
- Use a suitable lancet to prick the finger.
- Allow a good drop of blood to form before sampling.

Venous blood


- Patient samples stored refrigerated can be used without equilibration to room temperature.
- Mix the sample material well by inverting the tube 8-10 times before collecting a sample.

Afinion™ HbA1c Control

- Allow the control material to reach room temperature. (15-25°C) before use, which takes approximately 45 minutes.
- Mix the control material thoroughly by shaking the vial for 30 seconds.
- The sample can be extracted from the vial or the cap.

Filling the capillary

- Remove the sampling device from the Test Cartridge.
- Fill the capillary; bring the tip of the capillary just beneath the surface of the patient sample (a, b) or control material (c). Be sure that the capillary is completely filled, see arrow (d). It is not possible to overfill. Avoid air bubbles and excess sample on the outside of the capillary. Do not wipe off the capillary.
- Immediately replace the sampling device into the Test Cartridge.
- The analysis of the Test Cartridge should start within 2-3 minutes. If the Test Cartridge is stored too long before analysis, the sample material may dry or coagulate. An information code will then be displayed. Consult the Afinion™ AS100 Analyzer User Manual.
- Do not use a Test Cartridge that has been accidentally dropped on the floor or lab bench after sample collection.

ANALYSING A SAMPLE

- Analyse the Test Cartridge following the procedure described in the Afinion™ AS100 Analyzer User Manual.
- The analysis time is 3 minutes and 15 seconds.

INTERPRETATION OF RESULTS

Interpret the Afinion™ HbA1c test results with careful consideration to the patient's medical history, clinical examinations and other laboratory results. If the test result is questionable or if clinical signs and symptoms appear inconsistent with the test result, re-test the sample or confirm the result using another method. Analyse control materials frequently to verify the performance of the Afinion™ AS100 Analyzer System.

HbA1c measuring units

Three different measuring units are in use for reporting HbA1c test results¹¹:

- **mmol/mol** – the HbA1c values are aligned to the IFCC Reference Method for Measurements of HbA1c¹²
- **percentage (%)** – the HbA1c values are aligned to the assay used in the DCCT study, also known as NGSP-HbA1c^{3,13,14}
- **estimated average glucose (eAG)** – the HbA1c value is converted to an equivalent average glucose value¹⁵

Linear relationships have been established between the three units^{13,15}:

$$\text{NGSP-HbA1c (\%)} = \text{DCCT-HbA1c (\%)} = 0.092 \times \text{IFCC-HbA1c (mmol/mol)} + 2.15$$

$$\text{eAG (mmol/L)} = 1.59 \times \text{DCCT-HbA1c (\%)} - 2.59$$

Measuring range

The Afinion™ AS100 Analyzer displays the HbA1c value/concentration in mmol/mol, percentage (%), as the calculated estimated average glucose (eAG) or in combinations of these:

	HbA1c (%)	HbA1c mmol/mol	eAG mmol/L
Measuring range	4.0-15.0	20-140	3.8-21.3
Unit	0.1	1.0	0.1

 If the patient's HbA1c value is outside the measuring range, no test result will be reported and an information code will be displayed (see "Troubleshooting").

 **Important!** Patient's HbA1c values should be reported in units consistent with the national recommendations. Please contact your local supplier if the national recommendation is unknown. Consult the Afinion™ AS100 Analyzer User Manual for instructions on how to change the HbA1c measuring unit.

Standardisation

Afinion™ HbA1c is traceable to the IFCC Reference Method for Measurement of HbA1c⁵.

Reference range

For methods reporting DCCT traceable values, the upper limit of non-diabetic reference range is approximately 6% HbA1c^{7,8,9}.

Analytical specificity

Afinion™ HbA1c measures the total glycated hemoglobin and reports the HbA1c value.

The following hemoglobin (Hb) variants have been analysed and found not to affect the Afinion™ HbA1c test result: HbAC, HbAD, HbAE, HbF, HbAJ and HbAS¹⁰. Carbamylated and pre-glycated hemoglobin does not affect the Afinion™ HbA1c test result.

Interference

- No significant interference was observed at the following concentrations:

Bilirubin	342 µmol/L (20 mg/dL)
Triglycerides	15.7 mmol/L (1389 mg/dL)
Cholesterol	9.1 mmol/L (351 mg/dL)
Glucose	27.8 mmol/L (500 mg/dL)
Fructosamine	680 µmol/L

- | | |
|----------------------|------------------------|
| Acetaminophen | 1.7 mmol/L (256 µg/mL) |
| Ibuprofen | 1.8 mmol/L (372 µg/mL) |
| Acetylsalicylic acid | 3.3 mmol/L (599 µg/mL) |
| Salicylic acid | 4.3 mmol/L (593 µg/mL) |
| Glyburide | 3.9 µmol/L (256 µg/mL) |
| Metformin | 310 µmol/L |
- Hemolysis (in vitro) 5.0%
 - Anticoagulants (EDTA, heparin, citrate and NaF) at concentrations normally used in blood collection tubes do not interfere.

Important! It is possible that other substances and/or factors not listed above may interfere with the test and cause false results.

Limitations of the test

- Diluted samples cannot be used with Afinion™ HbA1c.
- Coagulated or hemolysed samples cannot be used with Afinion™ HbA1c.
- If the sample has a hemoglobin value below 6 g/dL or above 20 g/dL, no test result will be reported and an information code will be displayed.

QUALITY CONTROL

Quality control testing should be done to confirm that your Afinion™ AS100 Analyzer System is working properly and providing reliable results. Only when controls are used routinely and the values are within acceptable ranges can accurate results be assured for patient samples.

 It is recommended to keep a permanent record of all quality control results. The Afinion™ AS100 Analyzer automatically stores the control results in a separate log. Consult the Afinion™ AS100 Analyzer User Manual.

Choosing control material

 Afinion™ HbA1c Control from Axis-Shield is recommended for routine quality control testing. Consult the Afinion™ HbA1c Control Package Insert.

If controls from another supplier are used, the precision must be determined and acceptable ranges for the Afinion™ AS100 Analyzer System must be established. Lyophilised cell free lysate or hemolysed whole blood control materials cannot be used with Afinion™ HbA1c.

Frequency of control testing

Controls should be analysed:

- Anytime an unexpected test result is obtained.
- With each shipment of Afinion™ HbA1c test kits.
- With each new lot of Afinion™ HbA1c test kits.
- When training new operators in correct use of the Afinion™ HbA1c and the Afinion™ AS100 Analyzer.
- In compliance with national or local regulations.

Verifying the control results

 The measured value should be within the acceptable limits stated for the control material. Consult the Afinion™ HbA1c Control package insert.

If the result obtained for the Afinion™ HbA1c Control is outside the acceptable limits, make sure that:

- patient samples are not analysed until control results are within acceptable limits.
- the control vial has not passed its expiry date.
- the control vial has not been used for more than 60 days.
- the control vial and Afinion™ HbA1c Test Cartridge have been stored according to recommendations.
- there is no evidence of bacterial or fungal contamination of the control vial.

Correct any procedural error and re-test the control material.

If no procedural errors are detected:

- Examine the laboratory's quality control record to investigate the frequency of control failures.
- Ensure that there is no trend in out-of-range quality control results.
- Re-test the control material using a new control vial.
- Patient results must be declared invalid when controls do not perform as expected. Contact your local supplier for advice before analysing patient samples.

TROUBLESHOOTING

To ensure that correct HbA1c results are reported, the Afinion™ AS100 Analyzer performs optical, electronic and mechanical controls of the capillary, the Test Cartridge and all individual processing steps during the course of each analysis. When problems are detected by the built-in failsafe mechanisms, the Analyzer terminates the test and displays an information code.

The table below contains Afinion™ HbA1c specific information codes. Consult the Afinion™ AS100 Analyzer User Manual for information codes not listed in this table.

Code #	Cause
103	The hemoglobin concentration is below 6 g/dL
104	The hemoglobin concentration is above 20 g/dL
105	The HbA1c value is below measuring range
106	The HbA1c value is above measuring range

PERFORMANCE CHARACTERISTICS

Method comparison

A method comparison study, comprising 39 blood samples (4.9-11.7% HbA1c) were analysed with an affinity HPLC system and Afinion™ AS100 Analyzer. A linear regression analysis and a Bland-Altman analysis were performed. The results are shown in Table 1 and 2.

Table 1: Method comparison. Afinion™ HbA1c (y) vs. an affinity HPLC system (x). Linear regression analysis data. N=number of samples, r=correlation coefficient.

N	Regression line	r
39	$y = 0.96x + 0.33$	0.99

Table 2: Method comparison. Afinion™ HbA1c vs. an affinity HPLC system. Bland-Altman analysis. N=number of samples, CI=Confidence Interval.

N	Bias	95% CI
39	0.00% HbA1c	-0.31- 0.30% HbA1c

Precision

Within-run, between-day and total precision values were determined according to CLSI Protocol EP5-A. Afinion™ HbA1c Control C I, Control C II, one EDTA sample assayed for 20 days and one EDTA sample assayed for 17 days. The samples were analysed in duplicate twice a day. Precision data are summarised in Table 3.

Table 3: Within-run, between-day and total precision. N=number of days, CV=Coefficient of Variation.

Sample	N	Mean % HbA1c	Within-run CV (%)	Between-day CV (%)	Total CV (%)
Control C I	20	6.5	0.9	0.6	1.4
Control C II	20	9.1	0.6	0.5	1.1
Patient 1	17*	5.6	0.9	0.2	1.1
Patient 2	20	10.0	0.7	0.0	1.1

* Based on 17 days of analysis due to hemolysis of the sample.

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1. Lenzi S et al., The Clinical Usefulness of Glycated Hemoglobin in Monitoring Diabetes Mellitus: A long-Term Study. Clin Chem 1987; 33:55-56.
2. Frantzen F, Studies on glycated hemoglobin and new analytical methods for its determination. Dr. philos. thesis. NTNU, Trondheim 1998. ISBN 82-7861-116-5.
3. The Diabetes Control and Complications Trial Research Group, The Effect of Intensive Treatment of Diabetes on the Development and Progression of Long-Term Complications in Insulin-Dependent Diabetes Mellitus. N Engl J Med 1993; 329:977-986.
4. Sacks DB et al., Guidelines and Recommendations for Laboratory Analysis in the Diagnosis and Management of Diabetes Mellitus. Clin Chem 2002; 48:436-472.
5. Jeppsson, JO et al., Approved IFCC Reference Method for the Measurement of HbA1c in Human Blood, Clin Chem Lab Med 2002; 40(1):78-89.
6. Council Directive 67/548/EEC, Council Directive 76/769/EEC, Council Directive 1999/45/EC and Commission Directive 2001/60/EC.
7. Little, RR et al., Relationship of glycosylated hemoglobin to oral glucose tolerance. Implications for diabetes screening. Diabetes 1988;37(1):60-64.
8. Little, RR et al., Glycated haemoglobin predicts progression to diabetes mellitus in Pima Indians with impaired glucose tolerance. Diabetologia 1994; 37:252-256.
9. Rohlfing, CL et al., Use of GHb (HbA1c) in Screening for Undiagnosed Diabetes in the U.S. Population. Diabetes Care 2000; 23(2):187-191.
10. European Reference Laboratory, Manufacturer Check-Up Certification of Afinion™ HbA1c, 2004.
11. The American Diabetes Association, European Association for the Study of Diabetes, International Federation of Clinical Chemistry and Laboratory Medicine. Consensus Statement on the Worldwide Standardization of the Hemoglobin A1C Measurement. Diabetes Care 2007; 30:2399-2400.
12. Nordin G et al., Recommendation for term and measurement unit for "HbA1c". Clin Chem Lab Med 2007; 45:1081-1082.
13. Hoelzel W et al., IFCC reference system for measurement of haemoglobin A1c in human blood and the national standardization schemes in the United States, Japan, and Sweden: a method-comparison study. Clin Chem 2004; 50:166-174.
14. Weykamp C et al., The IFCC Reference Measurement System for HbA1c: a 6-year progress report. Clin Chem 2008; 54:24-248.
15. Nathan DM et al., A1c-Derived Average Glucose Study Group. Translating the A1C assay into estimated average glucose values. Diabetes Care 2008; 31:1473-1478

A distributor list is available at
www.axis-shield-poc.com

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