
REF 1115240

1115264 2007-05

A distributor list is available at [www.axis-shield-poc.com](http://www.axis-shield-poc.com)  
 Eine Distributorenliste kann unter [www.axis-shield-poc.com](http://www.axis-shield-poc.com) gefunden werden.  
 Un liste des distributeurs est disponible sur le site [www.axis-shield-poc.com](http://www.axis-shield-poc.com)  
 Per un elenco dei distributori, consultare [www.axis-shield-poc.com](http://www.axis-shield-poc.com)  
 En distributørliste finnes på [www.axis-shield-poc.com](http://www.axis-shield-poc.com)  
 En distributørlista finns på [www.axis-shield-poc.com](http://www.axis-shield-poc.com)  
 En distributørliste kan findes på [www.axis-shield-poc.com](http://www.axis-shield-poc.com)  
 Luettelo jakelijoista osoitteessa [www.axis-shield-poc.com](http://www.axis-shield-poc.com)

Target value and acceptable range/ Zielwert und toleranzbereich/ Valeur cible et gamme de mesure admissible/ Valore target ed intervallo di accettabilità/ Fasitverdi og akseptabelt område/ Facitvärde och acceptansgränser/ Målvärderi og acceptgrænser/ Tavoitearvo ja hyväksyttävä vaihteluväli

CONTROL C I

CONTROL C II

LOT

Albumin  
mg/L  
  
Creatinine  
mmol/L  
  
ACR  
mg/mmol

**Axis-Shield PoC AS**  
 P.O. Box 6863 Rodeløkka  
 NO-0504 Oslo, Norway  
[www.axis-shield-poc.com](http://www.axis-shield-poc.com)

ISO 9001 and ISO 13485 certified company

www.signatur.no

## SYMBOLS/ SYMBOLE/ SYMBOLES/ SIMBOLI/ SÍMBOLOS/ SYMBOLER/ SYMBOLIT

<p>Consult the user instructions          Bedienungsanleitungen beachten          Consultez les instructions d'utilisation          Consultare le istruzioni d'uso          Consulte las instrucciones para el usuario          Les bruksanvisningene          Läs bruksanvisningarna          Læs brugervejledningerne          Katso käyttöohjeet</p>	<p>Lot number          Chargen-Nr.          Numéro du lot          Numero di lotto          Número de lote          Lotnummer          Lotnummer          Lot-nummer          Eränumero</p>	<p>8°C 46°F  2°C 36°F</p> <p>Storage temperature 2-8°C (36-46°F)          Lagertemperatur 2-8°C (36-46°F)          Température de conservation 2-8°C (36-46°F)          Temperatura di conservazione 2-8°C (36-46°F)          Temperatura de almacenamiento 2-8°C (36-46°F)          Lagringstemperatur 2-8°C (36-46°F)          Förvaringstemperatur 2-8°C (36-46°F)          Opbevaringstemperatur 2-8°C (36-46°F)          Säilytyslämpötila 2-8°C (36-46°F)</p>
<p><i>In Vitro</i> Diagnostic Medical Device          Medizinisches Gerät zur <i>In-vitro</i>-Diagnostik          Appareil médical de diagnostic <i>in vitro</i>          Dispositivo medicale per la diagnostica <i>in vitro</i>          Dispositivo médico de diagnóstico <i>in vitro</i>          Medisinsk utstyr for <i>in vitro</i> diagnostikk  <i>In vitro</i> diagnostisk medicinsk utrustning  <i>In vitro</i> diagnostisk medicinsk udstyr  <i>In vitro</i> diagnostinen lääkinällinen laite</p>	<p>Biological risk          Biologisches Risiko          Risque biologique          Rischio biologico          Riesgo biológico          Biologisk risiko          Biologiskt riskmaterial          Biologisk risiko          Biologinen riski</p>	<p>Control C I          Kontrolle C I          Contrôle C I          Controllo C I          Control C I          Kontroll C I          Kontroll C I          Kontrol C I          Kontrolli C I</p>
<p>Catalogue number          Bestellnummer          Numéro de référence du catalogue          Numero di catalogo          Número de catálogo          Artikelnummer          Artikelnummer          Katalognummer          Luettelonumero</p>	<p>Manufacturer          Hersteller          Fabricant          Produttore          Fabricante          Produsent          Tillverkare          Producent          Valmistaja</p>	<p>Control C II          Kontrolle C II          Contrôle C II          Controllo C II          Control C II          Kontroll C II          Kontroll C II          Kontrol C II          Kontrolli C II</p>
<p>Expiry date (year-month)          Verfallsdatum (Jahr-Monat)          Date de péremption (année-mois)          Data di scadenza (anno-mese)          Fecha de vencimiento (año-mes)          Utløpsdato (år-måned)          Utgångsdatum (år-månad)          Udløbsdato (år-måned)          Viimeinen käyttöpäivämäärä (vuosi-kuukausi)</p>	<p>Conformity to the European directive 98/79/EC on <i>in vitro</i> diagnostic medical devices.          Konformität mit der europäischen Richtlinie 98/79/EC über medizinische Geräte zur <i>In-vitro</i>-Diagnostik.          Conformité avec la directive européenne 98/79/CE sur les appareils médicaux de diagnostic <i>in vitro</i>.          Conformità con la direttiva europea 98/79/CE sui dispositivi medici per la diagnostica <i>in vitro</i>.          Cumplimiento de la directiva europea 98/79/CE sobre los dispositivos médicos de diagnóstico <i>in vitro</i>.          Överensstemmelse med EU-direktiv 98/79/EC för <i>in vitro</i> diagnostisk medicinsk utstyr.          Överensstämmer med EU direktiv 98/79/EC för <i>in vitro</i> diagnostisk medicinsk utrustning.          Överensstemmelse med EU-direktivet 98/79/EC om <i>in vitro</i> diagnostisk udstyr til medicinsk brug.          Yhdenmukainen Euroopan direktiivin 98/79/EC kanssa koskien <i>in vitro</i> diagnostisia lääkinällisiä laitteita.</p>	

## BIBLIOGRAPHY/ LITERATUR/ BIBLIOGRAPHIE/ BIBLIOGRAFIA/ BIBLIOGRAFÍA/ LITTERATUR/ REFERENSER/ BIBLIOGRAFI/ KIRJALLISUUS

- ERM®, European Reference Material, [www.erm-crm.org](http://www.erm-crm.org).
- SRM, Standard Reference Material, National Institute for Standards and Technology, USA.

For use with Afinion™ ACR test / Afinion™ AS100 Analyzer

## PRODUCT DESCRIPTION

### Intended use

The Afinion™ ACR Control kit contains liquid preparations of albumin and creatinine in citrate buffer. The controls should be used to confirm that your Afinion™ AS100 Analyzer System is working properly and provides reliable results. Only when controls are used routinely and the values obtained are within acceptable ranges can accurate results for patient samples be assured.

### Kit contents

- Citrate buffer with albumin and creatinine at two concentration levels:
- 1 x 1,0 mL Afinion™ ACR Control C I
  - 1 x 1,0 mL Afinion™ ACR Control C II
  - 1 Package Insert

## WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use.
- Handle and dispose of the Afinion™ ACR Controls as potentially bio hazardous materials. Use gloves.
- The source material of Control C I and C II was found to be non-reactive for HbsAg and HIV.
- Do not use the Afinion™ ACR Controls after their expiry date or if they have not been stored in accordance with recommendations.
- Discard the vial if there is evidence of microbial or fungal contamination.

## STORAGE INSTRUCTIONS

- The expiry date of the kit only applies for unopened product stored refrigerated (2-8°C). The expiry date is the last day of the month stated on the outer container.
- Avoid exposure to direct sunlight and temperatures above 30°C.
- Do not freeze.

### Unopened control vials

- Unopened control vials are stable until expiry date stated on the vial label when stored refrigerated (2-8°C).

### Opened control vials

- Opened control vials are stable for 8 weeks when stored refrigerated (2-8°C). It is recommended to note the date of opening on the vial label.
- Always store the control vials refrigerated (2-8°C) in an upright position when not in use.
- The controls can be used directly from the refrigerator without equilibration to room temperature.

## FREQUENCY OF CONTROL TESTING

Good laboratory practice dictates that a quality control program should be established in all laboratories.

Controls should be analysed:

- anytime an unexpected test result is obtained.
- with each shipment of Afinion™ ACR test kits.
- with each new lot of Afinion™ ACR test kits.
- when training new operators in correct use of the Afinion™ ACR and the Afinion™ AS100 Analyzer.
- in compliance with national or local regulations.

 The Afinion™ AS100 Analyzer offers the possibility to store control test results in a separate log. Consult the Afinion™ AS100 User Manual for a detailed description.

## TESTING PROCEDURE

### Control sampling

- The controls can be used directly from the refrigerator without equilibration to room temperature.
- Mix the control material well by inverting the vial 8-10 times before collecting a sample.

 Consult the Afinion™ ACR Package Insert for instructions.

### Analysing a control

 Consult the Afinion™ AS100 Analyzer User Manual. The Afinion™ ACR Quick Guide also provides an illustrated step-by-step procedure.

### Verifying the control results

The Afinion™ ACR test used for target setting the Controls is calibrated against the reference solutions ERM®-DA470<sup>1</sup> (albumin) and SRM-914a<sup>2</sup> (creatinine). The measured values should be within the acceptable range stated for the Afinion™ ACR Control C I and C II on the front page.

If the result obtained is outside the acceptable range, make sure that:

- patient samples are not analysed until control results are within acceptable range.
- the control vial has not passed the expiry date.
- the control vial has not been used for more than 8 weeks.
- the control vial and the ACR test kit have been stored according to recommendations.
- there is no evidence of bacterial or fungal contamination of the control vial.

Correct any procedural error and re-test the control material.

If no procedural errors are detected:

- Re-test the control material using a new control vial.
- Examine the laboratory's quality control record to investigate the frequency of control failures.
- Ensure that there is no trend in out-of-range quality control results.
- Patient results must be declared invalid when controls do not perform as expected. Contact your local supplier for advice before analysing patient samples.

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