
Afinion™ HbA1c Control

REF 1115052

1115060 2004-12

A distributor list is available at www.axis-shield-poc.com
 Eine Distributorenliste kann unter www.axis-shield-poc.com gefunden werden.
 Un liste des distributeurs est disponible sur le site www.axis-shield-poc.com
 Per un elenco dei distributori, consultare il sito www.axis-shield-poc.com
 Una lista de distribuidores está disponible en www.axis-shield-poc.com
 En distributørliste finnes på www.axis-shield-poc.com
 En distributørlista finns på www.axis-shield-poc.com
 En distributørliste kan findes på www.axis-shield-poc.com
 Luettelo jakelijoista osoitteessa www.axis-shield.com

Expected values/ Erwartete Ergebnisse/ Valeurs attendues/
 Valori previsti/ Valores que se esperan/ Forventede verdier/
 Förväntade värden/ Forventede værdier/ Tavoitearvot

CONTROL C I

CONTROL C II

LOT

% HbA1c

% HbA1c SE

Axis-Shield PoC AS
 P.O. Box 6863 Rodeløkka
 NO-0504 Oslo, Norway
www.axis-shield-poc.com

ISO 9001 and ISO 13485 certified company

www.signalurno

SYMBOLS/ SYMBOLE/ SYMBOLES/ SIMBOLI/ SÍMBOLOS/ SYMBOLER/ SYMBOLIT

<p>Consult the Afinion™ user instructions Bedienungsanleitungen für Afinion™ beachten Consultez les instructions d'utilisation Afinion™ Consultare le istruzioni d'uso Afinion™ Consulte las instrucciones para el usuario de Afinion™ Les Afinion™ bruksanvisningene Läs bruksanvisningarna för Afinion™ Læs Afinion™ -brugervejledningerne Katso Afinion™ käyttöohjeet</p>	<p>Lot number Chargen-Nr. Numéro du lot Numero di lotto Número de lote Lotnummer Lotnummer Lot-nummer Eränumero</p>	<p>8°C 46°F 2°C 36°F</p> <p>Storage temperature 2-8°C (36-46°F) Lagertemperatur 2-8°C (36-46°F) Température de conservation 2-8°C (36-46°F) Temperatura di conservazione 2-8°C (36-46°F) Temperatura de almacenamiento 2-8°C (36-46°F) Lagringstemperatur 2-8°C (36-46°F) Förvaringstemperatur 2-8°C (36-46°F) Opbevaringstemperatur 2-8°C (36-46°F) Säilytys lämpötila 2-8°C (36-46°F)</p>
<p><i>In Vitro</i> Diagnostic Medical Device Medizinisches Gerät zur <i>In-vitro</i>-Diagnostik Appareil médical de diagnostic <i>in vitro</i> Dispositivo medicale per la diagnostica <i>in vitro</i> Dispositivo médico de diagnóstico <i>in vitro</i> Medisinsk utstyr til <i>in vitro</i> diagnostikk <i>In vitro</i> diagnostisk medicinsk utrustning <i>In vitro</i> diagnostisk medicinsk udstyr <i>In vitro</i> diagnostinen lääkinällinen laite</p>	<p>Biological risk Biologisches Risiko Risque biologique Rischio biologico Riesgo biológico Biologisk risiko Biologiskt riskmaterial Biologisk risiko Biologinen riski</p>	<p>Control C I Kontrolle C I Contrôle C I Controllo C I Control C I Kontroll C I Kontroll C I Control C I Kontrolli C I</p>
<p>Catalogue number Bestellnummer Numéro de référence du catalogue Numero di catalogo Número de catálogo Artikelnummer Artikelnummer Katalognummer Luettelunumero</p>	<p>Manufacturer Hersteller Fabricant Produttore Fabricante Produsent Tilverkare Producent Valmistaja</p>	<p>Control C II Kontrolle C II Contrôle C II Controllo C II Control C II Kontroll C II Kontroll C II Control C II Kontrolli C II</p>
<p>Expiration date (year-month) Verfallsdatum (Jahr-Monat) Date de péremption (année-mois) Data di scadenza (anno-mese) Fecha de vencimiento (año-mes) Utløpsdato (år-måned) Utgångsdatum (år-månad) Udløbsdato (år-måned) Viimeinen käyttöpäivämäärä (vuosi-kuukausi)</p>	<p>Conformity to the European directive 98/79/EC on <i>in vitro</i> diagnostic medical devices. Konformität mit der europäischen Richtlinie 98/79/EC über medizinische Geräte zur <i>In-vitro</i>-Diagnostik. Conformité avec la directive européenne 98/79/CE sur les appareils médicaux de diagnostic <i>in vitro</i>. Conformità con la direttiva europea 98/79/CE sui dispositivi medici per la diagnostica <i>in vitro</i>. Cumplimiento de la directiva europea 98/79/CE sobre los dispositivos médicos de diagnóstico <i>in vitro</i>. Överensstämelse med EU-direktiv 98/79/EC för <i>in vitro</i> diagnostisk medicinsk utstyr. Överensstämmer med EU direktiv 98/79/EC för <i>in vitro</i> diagnostisk medicinsk utrustning. Överensstemmelse med EU-direktiv 98/79/EC om <i>in vitro</i> diagnostisk udstyr til medicinsk brug. Yhdenmukainen Euroopan direktiivin 98/79/EC kanssa koskien <i>in vitro</i> diagnostisia lääkinällisiä laitteita.</p>	

BIBLIOGRAPHY/ LITERATUR/ BIBLIOGRAPHIE/ BIBLIOGRAFIA/ BIBLIOGRAFÍA/ LITTERATUR/ REFERENSER/ BIBLIOGRAFI/ KIRJALLISUUS

1. Council Directive 67/548/EEC, Council Directive 76/769/EEC, Council Directive 1999/45/EC and Commission Directive 2001/60/EC.
2. Jeppsson, JO et al., Approved IFCC Reference Method for the Measurement of HbA1c in Human Blood, Clin Chem Lab Med. 2002; 40(1):78-89.

For use with Afinion™ HbA1c / Afinion™ AS100 Analyzer.

PRODUCT DESCRIPTION

Intended use

Afinion™ HbA1c Control contains liquid preparations of stabilized porcine whole blood (Control C I) and human whole blood (Control C II). The Afinion™ HbA1c controls have been designed for use with the Afinion™ AS100 Analyzer System.

Quality control using the Afinion™ HbA1c Control should be done to confirm that your Afinion™ AS100 Analyzer System is working properly and provides reliable results. Only when controls are used routinely and the values obtained are within acceptable ranges can accurate results for patient samples be assured.

Kit contents

Control C I: Stabilized preparation from porcine whole blood (1 x 0,7 mL).
Control C II: Stabilized preparation from human whole blood (1 x 0,7 mL).

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use.
- Handle and dispose of the Afinion™ HbA1c controls as potentially bio hazardous materials. Use gloves.
- No known test method can offer complete assurance that products derived from human and animal blood will not transmit disease. Such materials should always be handled as potentially infectious.
- Control C II source material was found to be non-reactive for HbsAg, HCV, HIV-I and HIV-II when tested by FDA (Food and Drug Administration, USA) approved methods.
- The controls contain sodium azide as a preservative. The concentration is <0.1%, which is below that considered hazardous in normal use¹.
- Do not use the Afinion™ HbA1c controls after their expiry date or if they have not been stored in accordance with recommendations.
- Discard the vial if there is evidence of microbial or fungal contamination.

STORAGE INSTRUCTIONS

- The expiry date of the kit only applies if the product is stored at 2-8°C (36-46°F) in the original container. The expiry date is the last day of the month stated on the outer container.
- Avoid exposure to direct sunlight and temperatures above 30°C (86°F).
- Do not freeze.

Unopened control vials

- Unopened control vials are stable until expiry date indicated on the vial label when stored refrigerated (2-8°C, 36-46°F).

Opened control vials

- Opened control vials are stable for 30 days when stored refrigerated (2-8°C, 36-46°F). It is recommended to note the date of opening on the vial label.
- Always store the control vials refrigerated (2-8°C, 36-46°F) when not in use.
- Opened control vials should be stored in an upright position.

ANALYSING A CONTROL

Control sampling

Consult the Afinion™ HbA1c Package Insert for instructions.

Testing procedure

Consult the Afinion™ AS100 Analyzer User Manual for general Test Cartridge handling instructions and how to analyse a control on the Afinion™ AS100 Analyzer. An illustrated step-by-step procedure is also given in the Afinion™ HbA1c Quick Guide.

Interpretation of results

Good laboratory practice dictates that a quality control program should be established in all laboratories.

Consult the Afinion™ HbA1c Package Insert for recommended frequency of control testing for Afinion™ HbA1c.

The Afinion™ AS100 Analyzer offers the possibility to store control test results in a separate log. Consult the Afinion™ AS100 User Manual for a detailed description.

The measured values should be within the acceptable range for the control material. Re-test the control material if the control result is outside the range. Re-test a sample from a new control vial if the measured value is still outside the acceptable range. Patients' test results must be declared invalid if the control does not perform as expected. If the problem persists contact your local Afinion™ supplier for advice before analysing patient samples.

Expected values

Standardization of Afinion™ HbA1c is traceable to the IFCC Reference Method². The values reported are DCCT/NGSP aligned. In Sweden values traceable to Mono S are reported.

Two different values are stated in the package insert:

% HbA1c: Expected values according to the international standardization of Afinion™ HbA1c at DCCT/NGSP level, recommended in most countries.

% HbA1c SE: Expected values according to the Swedish Mono S reference method. For Swedish users only.

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