



Anticoagulant Clinic

Date :
 Number of Patients :
 Co-ordinator / Operator :

Clinic Set Up

Medical

- Mediswabs
- Sterile cotton wool balls
- Lancets
- Plasters
- Rubber Gloves
- Antiseptic spray (worktop cleaning)
- Clinical waste bin / sharps bin
- Syringes / needles / citrated sample tubes

Technical Equipment

- Thrombotrack™ Solo switch on
- Thrombotest™ reagent – prepare* (Batch No:)
- Calcium Chloride ()
- Normal QC Plasma – prepare* ()
- AK QC Plasma – prepare* ()
- Thrombotrack™ cuvettes, balls and ball dispenser
- Pipettes – set at 50ul, 250ul, 900ul and 4/3ml (for Thrombotest™ reagent)
- Pipette tips – yellow, blue and white
- Quick Reference Guide

* Reconstitute and prepare reagents as per Quick Reference Guide

Quality Control Results

Control Plasma Normal (INR range XXX-XXX)* Result – .

Control Plasma AK (INR rangeXXX -XXX)* Result – .

Quality control within range Yes/No

Comments

* insert lot specific range from product insert